

File 351:Derwent WPI 1963-2002/UD,UM &UP=200261

(c) 2002 Thomson Derwent

**\*File 351: Alerts can now have images sent via all delivery methods.**

See HELP ALERT and HELP PRINT for more info.

Set Items Description

-----

s pn=ru 2019189

S1 1 PN=RU 2019189

1/3,AB/1

DIALOG(R)File 351:Derwent WPI

(c) 2002 Thomson Derwent. All rts. reserv.

010258457

WPI Acc No: 1995-159712/ 199521

XRAM Acc No: C95-074164

**Type-specific Eimeria antigen prepn. - using an additional high-pressure liq. chromatography stage after the sterile oocysts have been homogenised and exposed to ultrasound**

Patent Assignee: LENG D EPIDEMIOLOGY MICROBIOL RES INST (LEEP-R)

Inventor: ARTYUKOV A I; MUNT YANU A I; VERBOV V N

Number of Countries: 001 Number of Patents: 001

Patent Family:

Patent No	Kind	Date	Applicat No	Kind	Date	Week
-----------	------	------	-------------	------	------	------

RU 2019189	C1	19940915	SU 5007122	A	19910802	199521 B
------------	----	----------	------------	---	----------	----------

Priority Applications (No Type Date): SU 5007122 A 19910802

Patent Details:

Patent No	Kind	Lan	Pg	Main IPC	Filing Notes
-----------	------	-----	----	----------	--------------

RU 2019189	C1	3	A61K-039/002		
------------	----	---	--------------	--	--

Abstract (Basic): RU 2019189 C

Eimeria antigens can be produced with increased activity and specificity by adding a high-pressure liq. chromatography stage to the process. After sterile sporulated oocysts have been obtd. in the usual way, then homogenised and exposed to ultrasound, sepn. is carried out using a grafted diol silica gel at 100 atmos., with an elution rate of 0.5 ml/min. at 4deg.C. The fractions with mol.wt. in the 100-200 kDa range are used as the antigen prepn.

USE - The method is used in veterinary science, esp. for obtaining antigen prepn.

ADVANTAGE - Crossover reactions with heterologous serums are reduced by a factor of 2-3 compared to the known technique.

Dwg.0/0